

## **EDTA Stock Solution (0.5M EDTA (pH 8.0) for 50 ml)**

Kathy Lu [sylu@ucsd.edu](mailto:sylu@ucsd.edu) 3/7/2013

Ref: [http://web.ku.edu/~sjmac/lab/reag/EDTA\\_0-5M\\_8-0\\_20061209.pdf](http://web.ku.edu/~sjmac/lab/reag/EDTA_0-5M_8-0_20061209.pdf)

1. 7.306 g EDTA [=  $C_{10}H_{14}O_8Na_2N_2 \cdot 2H_2O$  = EDTA disodium salt] + 40 ml dH<sub>2</sub>O
2. On shaker speed = 4 and mix
3. → Add 0.9-1g NaOH slowly until the solution is clear.
4. followed by 0.05-0.2ml 10M NaOH to pH8.0
5. → to 50 ml with dH<sub>2</sub>O